

Contextualized activity 1

Name: Maria Eduarda Franklin da Costa de Paula

Registration number: 2020020004

1.

Use of Lokomat

- Measure the patient's thigh segment;
- Get the waistcoat;
- If the waistcoat is aligned with User's Spine, then wear. Else, align it and wear them;
- Raise the patient on the platform;
- Lift the patient;
- Insert the patient in Lokomat's exoskeleton;
- Align the hip on exoskeleton;
- Align the legs on exoskeleton;
- Initialize the Lokomat;
- Set the time and Speed of the walk;
- While the time does not equal the 0, the patient walk;
- Show the process to the patient;
- If the time is equal 0 then, turn off the Lokomat;
- Decouple the exoskeleton;
- Put the patient on the wheelchair;
- Get the patient off the platform;
- Get off the patient's waistcoat.

Use of ZeroG

- Wear the waistcoat;
- Lower the support cable to the height of the User's head;
- If the waistcoat is dressed and the support cable is fitted, then lift the patient. Else, wear the waistcoat and fitted the support cable.
- If the Speed was adjusted the patient start walk, else, adjust them and next start walk.
- Repeat movement of ZeroG, while the user walk.

- If the person is unbalanced, the lift it;
- If the session time is over, then get patient off the ZeroG;
- Get off the waistcoat.

2. Lokomat

Memory:

Segment Size (Real), AlignmentWaistcoat (Boolean), DressedWaistcoat (Boolean), InsertedCorrectExoskeleton (Boolean), AlignmentHip (Boolean), AlignmentLegs (Boolean), LokomatInitialized (Boolean), time (Integer), Speed (Real), OnWheelchair (Boolean), InPlatform (Boolean), RaisedPatient (Boolean)

Processing:

MeasureSegment, AlignWaistcoat, GetWaistcoat, RaisePlatform, LeavePlatform, LiftPatient, InsertExoskeleton, AlignHip, AlignLegs, StartLokomat, SetTime, SetSpeed, StartWalk, ShowProcess, DecoupleExoskeleton, PutWheelchair, GetOffWaistcoat, wearWaistcoat;

Input:

Keyboard, Remote control, mouse, waistcoat's sensor.

Output:

Computer monitor, Led

ZeroG

Memory:

DressedWaistcoat (Boolean), FittedSupport-cable (Boolean), RaisedPatient (Boolean),

Speed Real), User walk (Boolean), Person-
Unbalanced (Boolean), SessionTime (Integer),
inZeroG (Boolean); SetSpeed (Boolean);

ZeroG

- Repeat movement of ZeroG, while the user walk.

Processing:

WearWaistcoat, PutSupportCable, Lift Pa-
tient, Set Speed, Patientwalk, ZeroG Move,
Set SessionTime, AnalyzeSessionTime,
GetOff ZeroG, GetOff Waistcoat, Put ZeroG.

Input:

Keyboard, mouse.

Output:

computer monitor, LED.

3. LoKomat

Structure of decision: LoKomat

- If the waistcoat is aligned with user's spine, then wear. Else, align and next wear them;
- If the time is equal 0 then, turn off the LoKomat.

ZeroG


- If the waistcoat is dressed and the Support is fitted, then lift the patient. Else wear them and fitted the cable;
- If the Speed was adjusted the patient start walk, else, adjust the speed and next start the walk;
- If the person is unbalanced, the lift it;
- If the session time is over, then get patient off the ZeroG.

Structure of repetition:



LoKomat:

- While the time does not equal the 0, the patient walk.

4. ZeroG



Lokomat



5.

Algorithm Automatic Lokomat

Var ClosedWaistcoat, FittedSupportCable,
ClosedExoskeleton, AlignmentHip,
AlignmentLegs Boolean

Var SessionTime, Integer

Var Speed, SizeSegment Real

Start

Input SizeSegment

If ClosedWaistcoat == True and FittedSupportCable == True then

Display "the cable will raise"

RaiseCable;

Else

If ClosedWaistcoat == False and FittedSupportCable == True then

display "the waistcoat is not closed"

Else if FittedCable == False and ClosedWaistcoat == True then

Display "the support cable is not fitted"

Else if FittedCable == False and ClosedWaistcoat == False then

Display "close the waistcoat and fit the support cable"

End-if

End-if

input ClosedExoskeleton, AlignmentHip, AlignmentLegs

If ClosedExoskeleton == True and AlignmentHip == True and AlignmentLegs == True then

 In: Input SessionTime, Speed

 While SessionTime != 0 do

 StartWalk

 SessionTime = SessionTime - 1;

 If SessionTime == 0 then

 Display "the session time is over, and support cable will down"

 DownCable;

 End-if

 End-While

Else

If ClosedExoskeleton == True and AlignmentHip == True and AlignmentLegs == False then
 Display "The Legs are not aligned"

Else if ClosedExoskeleton == True and AlignmentHip == False and AlignmentLegs == True then
 Display "the Hip is not Aligned"

Else if ClosedExoskeleton == False and AlignmentHip == True and AlignmentLegs == True then
 Display "The Lokomat's exoskeleton is not closed"

Else if ClosedExoskeleton == False and AlignmentHip == False and AlignmentLegs == False then
Display "It is necessary to align hip and legs. In addition, the exoskeleton
must be closed"

EndIf

End-if

End

Algorithm Automatic ZeroG

Var ClosedWaistcoat, FittedSupportCable, PersonUnbalanced, PatientWalk Boolean

Var SessionTime Integer

Var Speed Real

If ClosedWaistcoat == True and FittedSupportCable == True then

 input Speed, SessionTime

 While SessionTime != 0 do

 input PatientWalk

 if PatientWalk == False then

 StopWalk

 Else

 startWalk

 End-if

 SessionTime = SessionTime - 1

 if SessionTime == 0 then

 Display "the session time is over"

 StopWalk

 End-if

Else

 if ClosedWaistcoat == True and FittedSupportCable == False then

 Display "the support cable is not fitted"

 Else if ClosedWaistcoat == False and FittedSupportCable == True then

 Display "the waistcoat is not closed"

 Else if ClosedWaistcoat == False and FittedSupportCable == False then

 Display "Close the waistcoat and fit the support cable"

 End-if

End-if

End